



EIF2AK3 Antibody

Product Code	CSB-PA007512GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZJ5
Immunogen	Human PERK
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	eukaryotic translation initiation factor 2-alpha kinase 3;EIF2AK3;DKFZp781H1925;HRI;PEK;PERK;WRS ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF2AK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.