



EIF3G Antibody

Product Code	CSB-PA007536GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75821
Immunogen	Human EIF3G
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	eukaryotic translation initiation factor 3, subunit G;EIF3G;EIF3-P42;EIF3S4;eIF3-delta;eIF3-p44 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF3G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.