



EIF3I Antibody

Product Code	CSB-PA007538GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13347
Immunogen	Human EIF3I
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	eukaryotic translation initiation factor 3, subunit I;EIF3I;EIF3S2;PRO2242;TRIP-1;TRIP1;eIF3-beta;eIF3-p36 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF3I
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.