



EIF4E2 Antibody

Product Code	CSB-PA007558GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60573
Immunogen	Human EIF4E2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	eukaryotic translation initiation factor 4E family member 2;EIF4E2;4E-LP;4EHP;EIF4EL3;IF4e ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF4E2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.