



EIF4EBP2 Antibody

Product Code	CSB-PA007563GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13542
Immunogen	Human EIF4EBP2
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	eukaryotic translation initiation factor 4E binding protein 2; EIF4EBP2; 4EBP2; PHASII ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF4EBP2
Target Details	This gene encodes a member of the eukaryotic translation initiation factor 4E binding protein family. The gene products of this family bind eIF4E and inhibit translation initiation. However, insulin and other growth factors can release this inhibition via a phosphorylation-dependent disruption of their binding to eIF4E. Regulation of protein production through these gene products have been implicated in cell proliferation, cell differentiation and viral infection.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.