



EIF4G2 Antibody

Product Code	CSB-PA007569GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78344
Immunogen	Human EIF4G2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IP;WB:1:500-1:1000,IP:0.5-4.0ug
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	eukaryotic translation initiation factor 4 gamma, 2;EIF4G2;AAG1;DAP5;FLJ41344;NAT1;p97 ;
Immunogen Species	Homo sapiens (Human)
Target Names	EIF4G2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.