



ELK4 Antibody

Product Code	CSB-PA007607GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28324
Immunogen	Human ELK4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ELK4, ETS-domain protein (SRF accessory protein 1);ELK4;SAP1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ELK4
Target Details	This gene is a member of the Ets family of transcription factors and of the ternary complex factor (TCF) subfamily. Proteins of the TCF subfamily form a ternary complex by binding to the the serum response factor and the serum reponse element in the promoter of the c-fos proto-oncogene. This protein is phosphorylated by the kinases, MAPK1 and MAPK8. Several transcript variants have been described for this gene.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.