



ELOVL2 Antibody

Product Code	CSB-PA007620GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXB9
Immunogen	Human ELOVL2
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	elongation of very long chain fatty acids (FEN1/Elo2, SUR4/Elo3, yeast)-like 2; ELOVL2; FLJ20334; SSC2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ELOVL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.