



ENOPH1 Antibody

Product Code	CSB-PA007674GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHY7
Immunogen	Human ENOPH1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	enolase-phosphatase 1;ENOPH1;DKFZp586M0524;E1;FLJ12594;MASA;MST145 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ENOPH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.