



Recombinant Streptococcus suis Methionyl-tRNA formyltransferase (fmt)

Product Code	CSB-MP015150SUS
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4VTF3
Product Type	Recombinant Protein
Immunogen Species	Streptococcus suis (strain 05ZYH33)
Purity	≥85% (SDS-PAGE)
Sequence	MTKLIFMGTP AFSATVLRGL LADGNYNILA VVTQPDR AVG RKKVIQMTPV KEVALEYNLP VYQPEKLSGS QEMDELMNLG ADGIVTAAFG QFLPTKLLNS VDFAVNVHAS LLPKYRGGAP IHYALINGDE RAGVTIMEMV KEMDAGDMIS SDSIAIEESD NVGTLFEKLA VVGRDLLLQT LPAYIAGDLK PVAQNPEQVT FSPNIQPEEE VLDWNKTARQ LFNQIRGMYP WPVAHTYWQG ERFKIQEAVE AEGEGSVGRV IARSKKELII ATGQGALSLK TVQPAGKPKM TIADFLNGAG RDIAVGDQFG DQ
Source	Mammalian cell
Target Names	fmt
Protein Names	Recommended name: Methionyl-tRNA formyltransferase EC= 2.1.2.9
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.