



ENTPD8 Antibody

Product Code	CSB-PA007697GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5MY95
Immunogen	Human ENTPD8
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ectonucleoside triphosphate diphosphohydrolase 8; ENTPD8; GLSR2492; NTPDase-8; UNQ2492 ; NTPDase8
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ENTPD8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.