



EPCAM Antibody

Product Code	CSB-PA007717GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16422
Immunogen	Human EPCAM
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	epithelial cell adhesion molecule; EPCAM; CO-17A; CO17-1A; EGP; EGP-2; EGP34; EGP40; ESA; Ep-CAM; GA733-2; KS1/4; KSA; M4S1; MIC18; MK-1; TACST-1; TACSTD1; TROP1; hEGP-2 ; CD326
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EPCAM
Target Details	This gene encodes a carcinoma-associated antigen and is a member of a family that includes at least two type I membrane proteins. This antigen is expressed on most normal epithelial cells and gastrointestinal carcinomas and functions as a homotypic calcium-independent cell adhesion molecule. The antigen is being used as a target for immunotherapy treatment of human carcinomas. Mutations in this gene result in congenital tufting enteropathy.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.