



ERC2 Antibody

Product Code	CSB-PA007768GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15083
Immunogen	Human ERC2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ELKS/RAB6-interacting/CAST family member 2;ERC2;CAST;CAST1;ELKSL;KIAA0378;MGC133063;MGC133064;SPBC110; Spc110 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ERC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.