



ERCC6L Antibody

Product Code	CSB-PA007775GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2NKX8
Immunogen	Human ERCC6L
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	excision repair cross-complementing rodent repair deficiency, complementation group 6-like;ERCC6L;FLJ20105;MGC131695;PICH ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ERCC6L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.