



ERGIC3 Antibody

Product Code	CSB-PA007784GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y282
Immunogen	Human ERGIC3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ERGIC and golgi 3;ERGIC3;C20orf47;CGI-54;Erv46;NY-BR-84;PRO0989;SDBCAG84;dJ477O4.2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ERGIC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.