



ERO1L Antibody

Product Code	CSB-PA007797GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HE7
Immunogen	Human ERO1L
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ERO1-like (<i>S. cerevisiae</i>); ERO1L; ERO1-alpha ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ERO1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.