



ESRRG Antibody

Product Code	CSB-PA007837GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62508
Immunogen	Human ESRRG
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	estrogen-related receptor gamma;ESRRG;DKFZp781L1617;ERR3;FLJ16023;KIAA0832;NR3B3 ;ERRgamma
Product Type	Purified Rabbit Anti Human, Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ESRRG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.