



ETF1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA007840GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P62495 |
| Immunogen | Human ETF1 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | eukaryotic translation termination factor 1;ETF1;D5S1995;ERF;ERF1;MGC111066;RF1;SUP45L1;TB3-1 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | ETF1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |