



EVI2A Antibody

Product Code	CSB-PA007865GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22794
Immunogen	Human EVI2A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ecotropic viral integration site 2A;EVI2A;EVDA;EVI2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EVI2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.