



EXOSC3 Antibody

Product Code	CSB-PA007891GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQT5
Immunogen	Human EXOSC3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	exosome component 3;EXOSC3;CGI-102;MGC15120;MGC723;RP11-3J10.8;RRP40;Rrp40p;bA3J10.7;hRrp-40;hRrp40p;p10 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EXOSC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.