



EXTL3 Antibody

Product Code	CSB-PA007905GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43909
Immunogen	Human EXTL3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	exostoses (multiple)-like 3;EXTL3;DKFZp686C2342;EXTR1;KIAA0519;REG;REGR;RPR;botv ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EXTL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.