



Recombinant Human RNA polymerase II elongation factor ELL (ELL)

Product Code	CSB-EP007608HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P55199
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AALKEDRSY GLSCGRVSDG SKVSVFHVKL TDSALRAFES YRARQDSVSL RPSIRFQGSQ GHISIPQDC PAEARTFSFY LSNIGRDNPQ GSFDCIQQYV SSHGEVHLDC LGSIQDKITV CATDDSYQKA RQSMQAEEEE TRSRSAIVIK AGGRYLGKKV QFRKPAPGAT DAVPSRKRAT PINLASAIRK SGASAVSGGS GVSQRPFRDR VLHLLALRPY RKAELLRLQ KDGLTQADKD ALDGLLQQVA NMSAKDGTCT LQDCMYKDVQ KDWPGYSEGD QQLLKRVLVR KLCQPQSTGS LLGDPAASSP PGERGRSASP PQRKLQPPDF IDPLANKKPR ISHFTQRAQP AVNGKLGVPN GREALLPTPG PPASTDTLSS STHLPPRLEP PRAHDPLADV SNDLGHSGRD CEHGEEAAPA PTVRLGLPLL TDCAQPSRPH GSPSRSKPKK KSKKHKDKER AAEDKPRAQL PDCAPATHAT PGAPADTPGL NGTCSVSSVP TSTSETPDYL LKYAAISSSE QRQSYKNDFN AEYSEYRDLH ARIERITRRF TQLDAQLRQL SQGSEYETTT RGQILQEYRK IKKTNTNYSQ EKHRCEYLHS KLAHIKRLIA EYDQRQLQAW P
Source	E.coli
Target Names	ELL
Protein Names	Recommended name: RNA polymerase II elongation factor ELL Alternative name(s): Eleven-nineteen lysine-rich leukemia protein
Expression Region	2-621
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.