



# Recombinant Human RNA polymerase II elongation factor ELL (ELL)

<b>Product Code</b>	CSB-MP007608HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P55199
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AALKEDRSY GLSCGRVSDG SKVSVFHVKL TDSALRAFES YRARQDSVSL RPSIRFQGSQ GHISIPQPDC PAEARTFSFY LSNIGRDNPNQ GSFDCIQQYV SSHGEVHLDC LGSIQDKITV CATDDSYQKA RQSMQAEEEE TRSRSAIVIK AGGRYLGKKV QFRKPAGAT DAVPSRKRAT PINLASAIRK SGASAVSGGS GVSQRPFDR VLHLLALRPY RKAELLRLQ KDGLTQADKD ALDGLLQQVA NMSAKDGTCT LQDCMYKDVQ KDWPGYSEGD QLLKRVLVR KLCQPQSTGS LLGDPAASSP PGERGRSASP PQRKLPDF IDPLANKKPR ISHFTQRAQP AVNGKLGVPN GREALLPTPG PPASTDTLSS STHLPPRLEP PRAHDPLADV SNDLGHSGRD CEHGEEAAPA PTVRLGLPLL TDCAQPSRPH GSPSRSKPKK KSKKHDKER AAEDKPRAQL PDCAPATHAT PGAPADTPGL NGTCSVSSVP TSTSETPDYL LKYAAISSSE QRQSYKNDFN AEYSEYRDLH ARIERITRRF TQLDAQLRQL SQGSEYETT RGQILQEYRK IKKTNTNYSQ EKHRCEYLHS KLAHIKRLIA EYDQRQLQAW P
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ELL
<b>Protein Names</b>	Recommended name: RNA polymerase II elongation factor ELL Alternative name(s): Eleven-nineteen lysine-rich leukemia protein
<b>Expression Region</b>	2-621
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.