



FAF1 Antibody

Product Code	CSB-PA007963GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNN5
Immunogen	Human FAF1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	Fas (TNFRSF6) associated factor 1;FAF1;CGI-03;FLJ37524;HFAF1s;UBXD12;UBXN3A;hFAF1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FAF1
Target Details	Interaction of Fas ligand (TNFSF6) with the FAS antigen (TNFRSF6) mediates programmed cell death, also called apoptosis, in a number of organ systems. This protein binds to FAS antigen and can initiate apoptosis or enhance apoptosis initiated through FAS antigen. Initiation of apoptosis by This protein requires a ubiquitin-like domain but not the FAS-binding domain.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.