



FAM127A Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA008044GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O15255 |
| Immunogen | Human FAM127A |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA,IHC |
| Form | Liquid |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | RTL8C |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |