



FAM160B2 Antibody

Product Code	CSB-PA008108GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86V87
Immunogen	Human FAM160B2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	family with sequence similarity 160, member B2;FAM160B2;FLJ11125;FLJ21801;MGC138352;RAI16 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FAM160B2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.