



# FBXL18 Antibody

<b>Product Code</b>	CSB-PA008469GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96ME1
<b>Immunogen</b>	Human FBXL18
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,IHC
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	F-box and leucine-rich repeat protein 18;FBXL18;FLJ10776;FLJ11467;FLJ26934;FLJ38075;FLJ41541;Fbl18 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	FBXL18
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.