



# FEZ1 Antibody

<b>Product Code</b>	CSB-PA008600GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q99689
<b>Immunogen</b>	Human FEZ1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,IHC
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	fasciculation and elongation protein zeta 1 (zygin I);
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	FEZ1
<b>Target Details</b>	This gene is an ortholog of the <i>C. elegans</i> unc-76 gene, which is necessary for normal axonal bundling and elongation within axon bundles. Expression of this gene in <i>C. elegans</i> unc-76 mutants can restore to the mutants partial locomotion and axonal fasciculation, suggesting that it also functions in axonal outgrowth. The N-terminal half of the gene product is highly acidic. Alternatively spliced transcript variants encoding different isoforms of this protein have been described.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.