



FHL1 Antibody

Product Code	CSB-PA008663GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13642
Immunogen	Human Fhl1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	four and a half LIM domains 1; Fhl1; KyoT ;
Product Type	Purified Rabbit Anti Mus musculus PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FHL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.