



# Recombinant Mouse Alpha-2-macroglobulin receptor-associated protein (Lrpap1)

<b>Product Code</b>	CSB-EP013104MO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P55302
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QP IAGHGGKYSR EKNEPEMAAK RESGEEFRME KLNQLWEKAK RLHLSPVRLA ELHSDLKIQE RDELNWKCLK VEGLDKDGEK EAKLIHNLNV ILARYGLDGR KDAQMVHSNA LNEDTQDELG DPRLEKLWHK AKTSGKFSSE ELDKLWREFL HYKEKIQEYN VLLDTLSRAE EGYENLLSPS DMAHIKSDTL ISKHSELKDR LRSINQGLDR LRVSHQGYG STTEFEEPRV IDLWDLAQSA NFTEKELESF REELKHFEAK IEKHNHYQKQ LEISHQKCLKH VESIGDPEHI SRNKEYVLL EEKTKELGYK VKKHLQLDSS RVSRAHNL
<b>Source</b>	E.coli
<b>Target Names</b>	Lrpap1
<b>Protein Names</b>	Recommended name: Alpha-2-macroglobulin receptor-associated protein Short name= Alpha-2-MRAP Alternative name(s): Heparin-binding protein 44 Short name= HBP-44 Low density lipoprotein receptor-related protein-associated protein 1
<b>Expression Region</b>	29-360
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.