



FLT1 Antibody

Product Code	CSB-PA008732GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17948
Immunogen	Human FLT1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	fms-related tyrosine kinase 1 (vascular endothelial growth factor/vascular permeability factor receptor);FLT1;FLT;VEGFR1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FLT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.