



FOXG1 Antibody

Product Code	CSB-PA008815GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55316
Immunogen	Human FOXG1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	forkhead box G1;FOXG1;BF1;BF2;FHKL3;FKH2;FKHL1;FKHL2;FKHL3;FKHL4;FOXG1A;FOXG1B;FOXG1C;HBF-1;HBF-2;HBF-3;HBF- G2;HBF2;HFK1;HFK2;HFK3;KHL2;QIN ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FOXG1
Target Details	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in the development of the brain and telencephalon.