



# FOXG1 Antibody

<b>Product Code</b>	CSB-PA008815GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P55316
<b>Immunogen</b>	Human FOXG1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	forkhead box G1;FOXG1;BF1;BF2;FHKL3;FKH2;FKHL1;FKHL2;FKHL3;FKHL4;FOXG1A;FOXG1B;FOXG1C;HBF-1;HBF-2;HBF-3;HBF-G2;HBF2;HFK1;HFK2;HFK3;KHL2;QIN ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	FOXG1
<b>Target Details</b>	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in the development of the brain and telencephalon.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.