



Recombinant Horse Apolipoprotein A-II (APOA2)

| | |
|--------------------------|---|
| Product Code | CSB-EP001915HO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P83704 |
| Product Type | Recombinant Protein |
| Immunogen Species | Equus caballus (Horse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | QAEESCLQNL ASRYLQTVTD YGKDLVEKAL APELQAQAKA YFEKTQEQLT PLVKKIGNDL LNFFSHFIEL KTQPAT |
| Source | E.coli |
| Target Names | APOA2 |
| Protein Names | Recommended name: Apolipoprotein A-II Short name= Apo-AII Short name= ApoA-II Alternative name(s): Apolipoprotein A2 |
| Expression Region | 1-76 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Target Details | This gene encodes apolipoprotein (apo-) A-II, which is the second most abundant protein of the high density lipoprotein particles. The protein is found in plasma as a monomer, homodimer, or heterodimer with apolipoprotein D. Defects in this gene may result in apolipoprotein A-II deficiency or hypercholesterolemia. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |