



FRZB Antibody

Product Code	CSB-PA009008GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92765
Immunogen	Human FRZB
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	frizzled-related protein;FRZB;FRE;FRITZ;FRP-3;FRZB-1;FRZB-PEN;FRZB1;FZRB;SFRP3;SRFP3;hFIZ ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FRZB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.