



FXYD5 Antibody

Product Code	CSB-PA009093GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DB9
Immunogen	Human FXYD5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	FXYD domain containing ion transport regulator 5;FXYD5;HSPC113;IWU-1;IWU1;KCT1;OIT2;PRO6241;RIC;dysad ;; DYSADHERIN
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FXYD5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.