



FZR1 Antibody

Product Code	CSB-PA009113GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM11
Immunogen	Human FZR1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	fizzy/cell division cycle 20 related 1 (Drosophila);FZR1;CDC20C;CDH1;FZR;FZR2;HCDH;HCDH1;KIAA1242 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FZR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.