



GABARAPL1 Antibody

Product Code	CSB-PA009131GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0R8
Immunogen	Human GABARAPL1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	GABA(A) receptor-associated protein like 1;GABARAPL1;APG8-LIKE;APG8L;ATG8;ATG8L;GEC1 ;
Product Type	Purified Rabbit Anti human,mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GABARAPL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.