



GALNT15 Antibody

Product Code	CSB-PA009218GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N3T1
Immunogen	Human GALNTL2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GALNT15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.