



Recombinant Caldicellulosiruptor saccharolyticus Lipoyl synthase (lipA)

Product Code	CSB-EP012927CVX-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4XHV2
Product Type	Recombinant Protein
Immunogen Species	Caldicellulosiruptor saccharolyticus (strain ATCC 43494 / DSM 8903 / Tp8T 6331)
Purity	>85% (SDS-PAGE)
Sequence	MMQTKKPDWL RIRIKANQSV EEVIKLLKDL SLHTVCQEAQ CPNIFECFSK KTATFLILGD VCTRNCTFCD VKKGKQPQVD KNEPKKIAEA VKVLNLSYVV VTSVTRDDLE DGGAEHFANV IEKIKELNPQ TKVEVLIPDF NGDEKAIYKV VSARPDVLSH NVETVPRLYP TVRSKADYER SLSVLKVAKK MDNRIYTKSG LMVGLGETKE EVKEVLKNLR SVGCDFVTIG QYLSPSKQHY PVIEYIHPNV FEEYKEYAIS IGFKHVMSAP LVRSSYLAEE TTKII
Source	E.coli
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.