



GHITM Antibody

Product Code	CSB-PA009410GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3K2
Immunogen	Human GHITM
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	growth hormone inducible transmembrane protein;GHITM;DERP2;DKFZp566C0746;FLJ26584;HSPC282;MICS1;My021;P TD010;TMBIM5 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GHITM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.