



GIMAP5 Antibody

Product Code	CSB-PA009424GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96F15
Immunogen	Human GIMAP5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	GTPase, IMAP family member 5;GIMAP5;FLJ11296;HIMAP3;IAN4;IAN4L1;IAN5;IMAP3;hIAN5 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GIMAP5
Target Details	This gene encodes a protein belonging to the GTP-binding superfamily and to the immuno-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. In humans, the IAN subfamily genes are located in a cluster at 7q36.1.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.