



GJA3 Antibody

Product Code	CSB-PA009446GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6H8
Immunogen	Human GJA3
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	gap junction protein, alpha 3, 46kDa;GJA3;CX46;CZP3 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GJA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.