



GLIPR1L1 Antibody

Product Code	CSB-PA009504GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UWM5
Immunogen	Human GLIPR1L1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	GLI pathogenesis-related 1 like 1;GLIPR1L1;ALKN2972;MGC26856;PRO7434 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GLIPR1L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.