



GNAI2 Antibody

Product Code	CSB-PA009589GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04899
Immunogen	Human GNAI2
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC, IF, IP
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	guanine nucleotide binding protein (G protein), alpha inhibiting activity polypeptide 2;GNAI2;GIP;GNAI2B;H_LUCA15.1;H_LUCA16.1 ;
Immunogen Species	Homo sapiens (Human)
Target Names	GNAI2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.