



GNL2 Antibody

Product Code	CSB-PA009624GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13823
Immunogen	Human GNL2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	guanine nucleotide binding protein-like 2 (nucleolar);GNL2;FLJ40906;HUMAQUANTIG;NGP1;Ngp-1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GNL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.