



GORASP2 Antibody

Product Code	CSB-PA009676GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8Y8
Immunogen	Human GORASP2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	golgi reassembly stacking protein 2, 55kDa;GORASP2;DKFZp434D156;FLJ13139;GOLPH6;GRASP55;GRS2;p59 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GORASP2
Target Details	The Golgi complex plays a key role in the sorting and modification of proteins exported from the endoplasmic reticulum. This protein is involved in establishing the stacked structure of the Golgi apparatus.