



GPR142 Antibody

Product Code	CSB-PA009757GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z601
Immunogen	Human GPR142
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	G protein-coupled receptor 142;GPR142;PGR2 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GPR142
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.