



WLS Antibody

Product Code	CSB-PA009786GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T9L3
Immunogen	Human GPR177
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	G protein-coupled receptor 177;GPR177;C1orf139;EVI;FLJ23091;MGC131760;MGC14878;MRP;WLS ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	WLS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.