



GPR22 Antibody

Product Code	CSB-PA009796GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99680
Immunogen	Human GPR22
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	G protein-coupled receptor 22;GPR22;MGC129847 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GPR22
Target Details	This gene is a member of the G-protein coupled receptor 1 family and encodes a multi-pass membrane protein.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.